

For Research Use Only

# IRG1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68945-1-PBS



## Basic Information

Catalog Number:

68945-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28941

GenBank Accession Number:

NM\_001258406

GeneID (NCBI):

730249

UNIPROT ID:

A6NK06

Full Name:

immunoresponsive 1 homolog  
(mouse)

Calculated MW:

53 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2G12H4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50389-1, IRG1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68945-1-PBS. Detection antibody: 68945-2-PBS. Standard: Ag28941. Range: 0.781-100 ng/mL.