

For Research Use Only

IRG1 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 68945-2-PBS



Basic Information

Catalog Number: 68945-2-PBS	GenBank Accession Number: NM_001258406	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 730249	CloneNo.: 2G7E12
Source: Mouse	UNIPROT ID: A6NK06	
Isotype: IgG1	Full Name: immunoresponsive 1 homolog (mouse)	
Immunogen Catalog Number: AG28941	Calculated MW: 53 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50389-1, IRG1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68945-1-PBS. Detection antibody: 68945-2-PBS. Standard: Ag28941. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50389-2, IRG1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68945-2-PBS. Detection antibody: 68945-3-PBS. Standard: Ag28941. Range: 0.781-100 ng/mL.