

For Research Use Only

CREBL2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68949-1-PBS



Basic Information

Catalog Number:

68949-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG12492

GenBank Accession Number:

BC106052

GeneID (NCBI):

1389

UNIPROT ID:

O60519

Full Name:

cAMP responsive element binding
protein-like 2

Calculated MW:

120 aa, 14 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1C7B4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50393-1, CREBL2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68949-1-PBS. Detection antibody: 68949-2-PBS. Standard:Ag12492. Range: 0.391-100 ng/mL.



Cytometric bead array standard curve of MP50393-2, CREBL2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68949-1-PBS. Detection antibody: 68949-3-PBS. Standard:Ag12492. Range: 6.25-100 ng/mL.