

For Research Use Only

# CREBL2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68949-3-PBS



## Basic Information

Catalog Number:

68949-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG12492

GenBank Accession Number:

BC106052

GeneID (NCBI):

1389

UNIPROT ID:

O60519

Full Name:

cAMP responsive element binding  
protein-like 2

Calculated MW:

120 aa, 14 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3A4B9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50393-2, CREBL2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68949-1-PBS. Detection antibody: 68949-3-PBS. Standard: Ag12492. Range: 6.25-100 ng/mL