

For Research Use Only

# TBX1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68951-5-PBS



## Basic Information

Catalog Number: 68951-5-PBS	GenBank Accession Number: NM_005992	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 6899	CloneNo.: 2C11C6
Source: Mouse	UNIPROT ID: O43435	
Isotype: IgG2a	Full Name: T-box 1	
Immunogen Catalog Number: AG26792	Calculated MW: 43 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50397-3, TBX1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68951-5-PBS. Detection antibody: 68951-6-PBS. Standard: Ag26792. Range: 1.563-50 ng/mL