

For Research Use Only

PLIN1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68952-3-PBS



Basic Information

Catalog Number: 68952-3-PBS	GenBank Accession Number: NM_002666	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 5346	CloneNo.: 3F4D4
Source: Mouse	UNIPROT ID: O60240	
Isotype: IgG1	Full Name: perilipin	
Immunogen Catalog Number: AG34839	Calculated MW: 56 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50400-2, PLIN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68952-3-PBS. Detection antibody: 68952-2-PBS. Standard: Ag34839. Range: 1.563-50 ng/mL