

For Research Use Only

IL4I1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68955-2-PBS



Basic Information

Catalog Number: 68955-2-PBS	GenBank Accession Number: BC131625	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 259307	CloneNo.: 2E7E9
Source: Mouse	UNIPROT ID: Q96RQ9	
Isotype: IgG1	Full Name: interleukin 4 induced 1	
Immunogen Catalog Number: AG19136	Calculated MW: 567 aa, 63 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50424-1, Human IL41 Recombinant Matched Antibody Pair - PBS only. 68955-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag19136. 68955-2-PBS was HRP conjugated as the detection antibody. Range: 125-8000 pg/mL



Serum of sixteen individual healthy human donors was measured. The human IL41 concentration of detected samples was determined to be 629.0 pg/mL with a range of ND - 1913.6 pg/mL.