

For Research Use Only

# PAX8 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 68957-1-PBS



## Basic Information

Catalog Number:

68957-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG20348

GenBank Accession Number:

BC001060

GeneID (NCBI):

7849

UNIPROT ID:

Q06710

Full Name:

paired box 8

Calculated MW:

450 aa, 48 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1F11B1

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50408-1, PAX8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68957-1-PBS. Detection antibody: 68957-2-PBS. Standard:Ag20348. Range: 1.563-100 ng/mL.



Cytometric bead array standard curve of MP50408-2, PAX8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68957-3-PBS. Detection antibody: 68957-1-PBS. Standard:Ag20348. Range: 1.563-100 ng/mL.