

For Research Use Only

ISCU Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68958-3-PBS



Basic Information

Catalog Number:

68958-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6568

GenBank Accession Number:

BC061903

GeneID (NCBI):

23479

UNIPROT ID:

Q9H1K1

Full Name:

iron-sulfur cluster scaffold homolog
(E. coli)

Calculated MW:

18 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2G9H5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50409-2, ISCU Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68958-3-PBS. Detection antibody: 68958-4-PBS. Standard: Ag6568. Range: 0.781-100 ng/mL.