

For Research Use Only

MBTPS2 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 68959-2-PBS



Basic Information

Catalog Number: 68959-2-PBS	GenBank Accession Number: BC036465	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 51360	CloneNo.: 3D1A5
Source: Mouse	UNIPROT ID: O43462	
Isotype: IgG2a	Full Name: membrane-bound transcription factor peptidase, site 2	
Immunogen Catalog Number: AG4390	Calculated MW: 36 kDa, 57 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50410-1, MBTPS2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68959-1-PBS. Detection antibody: 68959-2-PBS. Standard:Ag4390. Range: 0.195-100 ng/mL.



Cytometric bead array standard curve of MP50410-3, MBTPS2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68959-2-PBS. Detection antibody: 68959-3-PBS. Standard:Ag4390. Range: 0.391-100 ng/mL.