

For Research Use Only

COX6A1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68960-1-PBS



Basic Information

Catalog Number:

68960-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG2006

GenBank Accession Number:

BC007723

GeneID (NCBI):

1337

UNIPROT ID:

P12074

Full Name:

cytochrome c oxidase subunit VIa

polypeptide 1

Calculated MW:

85 aa, 9 kDa

Purification Method:

Protein G purification

CloneNo.:

1F5D5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50413-1, COX6A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68960-1-PBS. Detection antibody: 68960-2-PBS. Standard: Ag2006. Range: 0.098-25 ng/mL



Cytometric bead array standard curve of MP50413-2, COX6A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68960-1-PBS. Detection antibody: 68960-3-PBS. Standard: Ag2006. Range: 0.098-25 ng/mL