

For Research Use Only

# COX6A1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68960-2-PBS



## Basic Information

Catalog Number:

68960-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG2006

GenBank Accession Number:

BC007723

GeneID (NCBI):

1337

UNIPROT ID:

P12074

Full Name:

cytochrome c oxidase subunit VIa polypeptide 1

Calculated MW:

85 aa, 9 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2E7E9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50413-1, COX6A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68960-1-PBS. Detection antibody: 68960-2-PBS. Standard:Ag2006. Range: 0.098-25 ng/mL



Cytometric bead array standard curve of MP50413-4, COX6A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68960-5-PBS. Detection antibody: 68960-2-PBS. Standard:Ag2006. Range: 0.098-25 ng/mL