

For Research Use Only

COX6A1 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 68960-5-PBS



Basic Information

Catalog Number: 68960-5-PBS	GenBank Accession Number: BC007723	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 1337	CloneNo.: 1G11H2
Source: Mouse	UNIPROT ID: P12074	
Isotype: IgG1	Full Name: cytochrome c oxidase subunit VIa polypeptide 1	
Immunogen Catalog Number: AG2006	Calculated MW: 85 aa, 9 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50413-3, COX6A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68960-4-PBS. Detection antibody: 68960-5-PBS. Standard:Ag2006. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50413-4, COX6A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68960-5-PBS. Detection antibody: 68960-2-PBS. Standard:Ag2006. Range: 0.098-25 ng/mL.