

For Research Use Only

# GPR64 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68963-5-PBS



## Basic Information

Catalog Number:

68963-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG31274

GenBank Accession Number:

BC099901

GeneID (NCBI):

10149

UNIPROT ID:

Q8IZP9

Full Name:

G protein-coupled receptor 64

Calculated MW:

1017 aa, 112 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1B11H3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50415-4, GPR64 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68963-5-PBS. Detection antibody: 68963-6-PBS. Standard: Ag31274. Range: 3.125-100 ng/mL.