

For Research Use Only

# p38 MAPK Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68965-1-PBS



## Basic Information

Catalog Number:

68965-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC031574

GeneID (NCBI):

1432

UNIPROT ID:

Q16539

Full Name:

mitogen-activated protein kinase 14

Calculated MW:

360 aa, 41 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1A11D8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50419-1, p38 MAPK Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66234-1-PBS. Detection antibody: 68965-1-PBS. Standard: Ag5797. Range: 0.195-200 ng/mL