

For Research Use Only

# p38 MAPK phospho-T180/Y182 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68966-1-PBS



## Basic Information

Catalog Number:

68966-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

NM\_139014

GeneID (NCBI):

1432

UNIPROT ID:

Q16539

Full Name:

mitogen-activated protein kinase 14

Purification Method:

Protein G Magarose purification

CloneNo.:

2B5D5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50419-2, p38 MAPK Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66234-1-PBS. Detection antibody: 68966-1-PBS. Standard: Ag5797. Range: 1.563-200 ng/mL