

For Research Use Only

Integrin alpha-6 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 68968-2-PBS



Basic Information

Catalog Number: 68968-2-PBS	GenBank Accession Number: BC136456	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 3655	CloneNo.: 2A9B3
Source: Mouse	UNIPROT ID: P23229	
Isotype: IgG2a	Full Name: integrin, alpha 6	
Immunogen Catalog Number: AG16335	Calculated MW: 1073 aa, 119 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50396-4, Integrin alpha-6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68968-1-PBS. Detection antibody: 68968-2-PBS. Standard: Ag16335. Range: 0.195-100 ng/mL.



Cytometric bead array standard curve of MP50396-5, Integrin alpha-6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68968-2-PBS. Detection antibody: 68968-3-PBS. Standard: Ag16335. Range: 0.781-100 ng/mL.