

For Research Use Only

NDUFS4 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68971-2-PBS



Basic Information

Catalog Number: 68971-2-PBS	GenBank Accession Number: BC005270	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 4724	CloneNo.: 1G11D11
Source: Mouse	UNIPROT ID: O43181	
Isotype: IgG1	Full Name: NADH dehydrogenase (ubiquinone) Fe-S protein 4, 18kDa (NADH- coenzyme Q reductase)	
Immunogen Catalog Number: AG8566	Calculated MW: 175 aa, 20 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50426-1, NDUF54 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68971-1-PBS. Detection antibody: 68971-2-PBS. Standard: Ag8566. Range: 0.391-50 ng/mL