

For Research Use Only

# IMMP2L Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 68972-2-PBS



## Basic Information

Catalog Number:

68972-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG8752

GenBank Accession Number:

BC008497

GeneID (NCBI):

83943

UNIPROT ID:

Q96T52

Full Name:

IMP2 inner mitochondrial membrane  
peptidase-like (S. cerevisiae)

Calculated MW:

12 kDa, 20 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1F4B4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50430-1, IMMP2L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68972-1-PBS. Detection antibody: 68972-2-PBS. Standard: Ag8752. Range: 0.195-3.125 ng/mL.



Cytometric bead array standard curve of MP50430-2, IMMP2L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68972-2-PBS. Detection antibody: 68972-3-PBS. Standard: Ag8752. Range: 0.098-50 ng/mL.