

For Research Use Only

IMMP2L Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68972-3-PBS



Basic Information

Catalog Number:

68972-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG8752

GenBank Accession Number:

BC008497

GeneID (NCBI):

83943

UNIPROT ID:

Q96T52

Full Name:

IMP2 inner mitochondrial membrane
peptidase-like (S. cerevisiae)

Calculated MW:

12 kDa, 20 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1H7D6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50430-2, IMMP2L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68972-2-PBS. Detection antibody: 68972-3-PBS. Standard: Ag8752. Range: 0.098-50 ng/mL