

For Research Use Only

CMPK2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68974-1-PBS



Basic Information

Catalog Number:

68974-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG22885

GenBank Accession Number:

BC141802

GeneID (NCBI):

129607

UNIPROT ID:

Q5EBM0

Full Name:

cytidine monophosphate (UMP-CMP)
kinase 2, mitochondrial

Calculated MW:

449 aa, 49 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2H5C11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50436-1, CMPK2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68974-1-PBS. Detection antibody: 68974-2-PBS. Standard: Ag22885. Range: 0.391-50 ng/mL



Cytometric bead array standard curve of MP50436-2, CMPK2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68974-1-PBS. Detection antibody: 68974-3-PBS. Standard: Ag22885. Range: 0.391-50 ng/mL