

For Research Use Only

HMGCR Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68975-3-PBS



Basic Information

Catalog Number:

68975-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG4444

GenBank Accession Number:

BC033692

GeneID (NCBI):

3156

UNIPROT ID:

P04035

Full Name:

3-hydroxy-3-methylglutaryl-
Coenzyme A reductase

Calculated MW:

888 aa, 95 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3H3E3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50438-2, HMGR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68976-1-PBS. Detection antibody: 68975-3-PBS. Standard: Ag4444. Range: 0.098-25 ng/mL