

For Research Use Only

OPN, SPP1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68983-2-PBS



Basic Information

Catalog Number: 68983-2-PBS	GenBank Accession Number: BC007016	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 6696	CloneNo.: 1E12B12
Source: Mouse	UNIPROT ID: P10451-5	
Isotype: IgG1	Full Name: secreted phosphoprotein 1	
Immunogen Catalog Number: AG19216	Calculated MW: 314 aa, 35 kDa	

Applications

Tested Applications:
IHC, Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50455-1, OPN, SPP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68983-1-PBS. Detection antibody: 68983-2-PBS. Standard: Ag19216. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50455-4, OPN, SPP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68983-4-PBS. Detection antibody: 68983-2-PBS. Standard: Ag19216. Range: 3.125-100 ng/mL.