

For Research Use Only

PTGES Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68984-1-PBS



Basic Information

Catalog Number: 68984-1-PBS	GenBank Accession Number: BC008280	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 9536	CloneNo.: 1E7F1
Source: Mouse	UNIPROT ID: O14684	
Isotype: IgG1	Full Name: prostaglandin E synthase	
Immunogen Catalog Number: AG33556	Calculated MW: 17 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP50456-1, PTGES Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68984-1-PBS. Detection antibody: 68984-2-PBS. Standard:Ag33556. Range: 0.195-100 ng/mL



Cytometric bead array standard curve of MP50456-2, PTGES Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68984-1-PBS. Detection antibody: 68984-3-PBS. Standard:Ag33556. Range: 0.195-100 ng/mL



Cytometric bead array standard curve of MP50456-3, PTGES Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68984-1-PBS. Detection antibody: 68984-4-PBS. Standard:Ag33556. Range: 0.195-100 ng/mL