

For Research Use Only

Connexin-26 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68991-1-PBS



Basic Information

Catalog Number: 68991-1-PBS	GenBank Accession Number: BC017048	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 2706	CloneNo.: 1F8G5
Source: Mouse	UNIPROT ID: P29033	
Isotype: IgG1	Full Name: gap junction protein, beta 2, 26kDa	
Immunogen Catalog Number: AG10666	Calculated MW: 26 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50480-1, Connexin-26 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68991-1-PBS. Detection antibody: 68991-2-PBS. Standard: Ag10666. Range: 1.563-100 ng/mL.