

For Research Use Only

LAPTM4A Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68994-1-PBS



Basic Information

Catalog Number:

68994-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG6768

GenBank Accession Number:

BC000421

GeneID (NCBI):

9741

UNIPROT ID:

Q15012

Full Name:

lysosomal protein transmembrane 4

alpha

Calculated MW:

27 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2E5B3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50487-1, LAPTM4A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68994-1-PBS. Detection antibody: 68994-2-PBS. Standard: Ag6768. Range: 6.25-100 ng/mL.