

For Research Use Only

NUP210 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68996-1-PBS



Basic Information

Catalog Number: 68996-1-PBS	GenBank Accession Number: BC067089	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 23225	CloneNo.: 1C6A1
Source: Mouse	UNIPROT ID: Q8TEM1	
Isotype: IgG2a	Full Name: nucleoporin 210kDa	
Immunogen Catalog Number: AG10494	Calculated MW: 967aa,106 kDa; 1887aa,205 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50488-1, NUP210 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68996-1-PBS. Detection antibody: 68996-2-PBS. Standard:Ag10494. Range: 0.391-100 ng/mL.