

For Research Use Only

# IL-10 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 69018-2-PBS

1 Publications

## Basic Information

Catalog Number:

69018-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

HZ-1145

GenBank Accession Number:

GeneID (NCBI):

3586

Full Name:

interleukin 10

Purification Method:

Protein A purification

CloneNo.:

7B8G3

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Cited Applications:

IHC

Species Specificity:

Human

Cited Species:

rat

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP50035-1, IL-10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 69018-1-PBS. Detection antibody: 69018-2-PBS. Standard: HZ-1145. Range: 15.6-1000 pg/mL.