

For Research Use Only

mouse IgG Kappa chain Recombinant antibody, PBS Only



Catalog Number:80015-1-PBS

Basic Information

Catalog Number:

80015-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Calculated MW:

155 kDa

Purification Method:

Protein A purification

CloneNo.:

3D9

Applications

Tested Applications:

WB,Indirect ELISA

Species Specificity:

human, mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HEK-293 cells were subjected to SDS-PAGE followed by western blot with eight different monoclonal antibodies and subsequently detected by recombinant rabbit anti-mouse kappa chain antibody 80015-1-Ig at 1:10000 and then probed with HRP-goat anti-rabbit secondary antibody. This data was developed using the same antibody clone with 80015-1-PBS in a different storage buffer formulation.