

For Research Use Only

STAT2 Recombinant antibody, PBS Only (Detector)

Catalog Number:80061-4-PBS



Basic Information

Catalog Number:

80061-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10168

GenBank Accession Number:

BC051284

GeneID (NCBI):

6773

UNIPROT ID:

P52630

Full Name:

signal transducer and activator of
transcription 2, 113kDa

Calculated MW:

851 aa, 98 kDa

Purification Method:

Protein A purification

CloneNo.:

242685E6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01853-1, STAT2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80061-3-PBS. Detection antibody: 80061-4-PBS. Standard: Ag10168. Range: 1.563-200 ng/mL.



Cytometric bead array standard curve of MP01853-2, STAT3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80061-2-PBS. Detection antibody: 80061-4-PBS. Standard: Ag10168. Range: 1.563-200 ng/mL.