

For Research Use Only

TAU Recombinant antibody, PBS Only (Capture)



Catalog Number: 80073-3-PBS

Basic Information

Catalog Number:

80073-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC000558

GeneID (NCBI):

4137

UNIPROT ID:

P10636

Full Name:

microtubule-associated protein tau

Calculated MW:

37-46, 79-81 kDa

Purification Method:

Protein A purification

CloneNo.:

230429G9

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00210-3, TAU Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80073-3-PBS. Detection antibody: 80073-5-PBS. Standard: SY01799P. Range: 0.156-2.5 ng/mL