

For Research Use Only

Mouse Tnfa Recombinant antibody, PBS Only (Detector)



Catalog Number: 80258-2-PBS

Basic Information

Catalog Number:
80258-2-PBS

Size:
1 mg/ml

Source:
Rabbit

Isotype:
IgG

GenBank Accession Number:
NM-013693

GeneID (NCBI):
21926

UNIPROT ID:
P06804

Full Name:
tumor necrosis factor

Purification Method:
Protein A purification

CloneNo.:
240702B3

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Mouse

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00682-1, MOUSE TNF-alpha Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80258-3-PBS. Detection antibody: 80258-2-PBS. Standard: Eg31366. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00682-2, MOUSE TNF-alpha Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80258-4-PBS. Detection antibody: 80258-2-PBS. Standard: Eg31366. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP00682-3, MOUSE TNF-alpha Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80258-5-PBS. Detection antibody: 80258-2-PBS. Standard: Eg31366. Range: 0.313-40 ng/mL