

For Research Use Only



Mouse TNF-alpha Recombinant antibody, PBS Only (Capture)

Catalog Number:80258-4-PBS

Basic Information

Catalog Number:

80258-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM-013693

GeneID (NCBI):

21926

UNIPROT ID:

P06804

Full Name:

tumor necrosis factor

Purification Method:

Protein A purification

CloneNo.:

240702C9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00682-2, MOUSE TNF-alpha Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80258-4-PBS. Detection antibody: 80258-2-PBS. Standard: Eg31366. Range: 0.156-20 ng/mL.