

For Research Use Only

XRCC5 Recombinant antibody, PBS Only (Detector)



Catalog Number: 80325-4-PBS

Basic Information

Catalog Number: 80325-4-PBS	GenBank Accession Number: BC019027	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 7520	CloneNo.: 240301B10
Source: Rabbit	UNIPROT ID: P13010	
Isotype: IgG	Full Name: X-ray repair complementing defective repair in Chinese hamster cells 5 (double-strand-break rejoining)	
Immunogen Catalog Number: AG9454	Calculated MW: 732 aa, 83 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00384-1, XRCC5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80325-5-PBS. Detection antibody: 80325-4-PBS. Standard: Ag9454. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00384-3, XRCC5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80325-3-PBS. Detection antibody: 80325-4-PBS. Standard: Ag9454. Range: 0.625-80 ng/mL