

For Research Use Only

XRCC5 Recombinant antibody, PBS Only



Catalog Number: 80325-6-PBS

Basic Information

Catalog Number: 80325-6-PBS	GenBank Accession Number: BC019027	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 7520	CloneNo.: 240301G6
Source: Rabbit	UNIPROT ID: P13010	
Isotype: IgG	Full Name: X-ray repair complementing defective repair in Chinese hamster cells 5 (double-strand-break rejoining)	
Immunogen Catalog Number: AG9454	Calculated MW: 732 aa, 83 kDa	

Applications

Tested Applications:
WB, IF, ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 80325-6-RR (XRCC5 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 80325-6-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using XRCC5 antibody (80325-6-RR, Clone: 240301G6) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 80325-6-PBS in a different storage buffer formulation.