

For Research Use Only

# ADRP Recombinant antibody, PBS Only (Detector)



Catalog Number: 80362-4-PBS

## Basic Information

Catalog Number:

80362-4-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7539

GenBank Accession Number:

BC005127

GeneID (NCBI):

123

UNIPROT ID:

Q99541

Full Name:

adipose differentiation-related protein

Calculated MW:

48 kDa

Purification Method:

Protein A purification

CloneNo.:

230274E5

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP00071-1, ADRP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80362-3-PBS. Detection antibody: 80362-4-PBS. Standard: Ag7539. Range: 0.156-10 ng/mL