

For Research Use Only

E-cadherin Recombinant antibody, PBS Only (Detector)



Catalog Number: 80541-6-PBS

Basic Information

Catalog Number:

80541-6-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14973

GenBank Accession Number:

BC141838

GeneID (NCBI):

999

UNIPROT ID:

P12830

Full Name:

cadherin 1, type 1, E-cadherin
(epithelial)

Calculated MW:

882 aa, 97 kDa

Purification Method:

Protein A purification

CloneNo.:

240159E2

Applications

Tested Applications:

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP00187-2, Human E-cadherin Recombinant Matched Antibody Pair - BSA and Azide Free measured by Cytometric bead array. Capture antibody: 80541-4-PBS. Detection antibody: 80541-6-PBS. Standard: Ag14973. Range: 0.625-20 ng/mL



Standard curve of MP00187-3, Human E-cadherin Recombinant Matched Antibody Pair - BSA and Azide Free measured by Cytometric bead array. Capture antibody: 80541-3-PBS. Detection antibody: 80541-6-PBS. Standard: Ag14973. Range: 1.25-20 ng/mL