

For Research Use Only

Tubulin-beta Recombinant antibody, PBS Only (Detector)



Catalog Number:80676-1-PBS

Basic Information

Catalog Number: 80676-1-PBS	GenBank Accession Number: BC000748	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 10381	CloneNo.: 8P15
Source: Rabbit	UNIPROT ID: Q13509	
Isotype: IgG	Full Name: tubulin, beta 3	
Immunogen Catalog Number: AG0117	Calculated MW: 450 aa, 50 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00550-1, Tubulin-beta Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80680-1-PBS. Detection antibody: 80676-1-PBS. Standard: ag0117. Range: 0.098-100 ng/mL.