

For Research Use Only

PTEN Recombinant antibody, PBS Only (Capture)



Catalog Number: 80718-5-PBS

Basic Information

Catalog Number: 80718-5-PBS	GenBank Accession Number: BC005821	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 5728	CloneNo.: 240134B11
Source: Rabbit	UNIPROT ID: P60484	
Isotype: IgG	Full Name: phosphatase and tensin homolog	
Immunogen Catalog Number: AG17274	Calculated MW: 47 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00328-1, PTEN Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80718-5-PBS. Detection antibody: 80718-4-PBS. Standard: Ag17274. Range: 0.313-20 ng/mL



Cytometric bead array standard curve of MP00328-2, PTEN Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80718-5-PBS. Detection antibody: 80718-3-PBS. Standard: Ag17274. Range: 0.313-20 ng/mL