

For Research Use Only

Tubulin-Alpha Recombinant antibody



Catalog Number: 80740-1-RR

Basic Information

Catalog Number: 80740-1-RR	GenBank Accession Number: BC009314	Purification Method: Protein A purification
Size:	GeneID (NCBI): 10376	CloneNo.: 3H19
Source: Rabbit	UNIPROT ID: P68363	
Isotype: IgG	Full Name: tubulin, alpha 1b	
Immunogen Catalog Number: AG18034	Calculated MW: 50 kDa	

Applications

Tested Applications:
ELISA
Species Specificity:
Human

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data