

For Research Use Only

p38 MAPK Recombinant antibody, PBS Only (Capture)



Catalog Number: 80821-5-PBS

Basic Information

Catalog Number: 80821-5-PBS	GenBank Accession Number: BC031574	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1432	CloneNo.: 240181E7
Source: Rabbit	UNIPROT ID: Q16539	
Isotype: IgG	Full Name: mitogen-activated protein kinase 14	
Immunogen Catalog Number: AG5115	Calculated MW: 360 aa, 41 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00437-1, p38 MAPK Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80821-5-PBS. Detection antibody: 80821-6-PBS. Standard: Ag5115. Range: 1.56-100 ng/mL



Cytometric bead array standard curve of MP00437-2, p38 MAPK Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80821-5-PBS. Detection antibody: 80821-4-PBS. Standard: Ag5115. Range: 3.13-100 ng/mL