

For Research Use Only

# p38 MAPK Recombinant antibody, PBS Only (Detector)



Catalog Number:80821-6-PBS

## Basic Information

Catalog Number: 80821-6-PBS	GenBank Accession Number: BC031574	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1432	CloneNo.: 240181H8
Source: Rabbit	UNIPROT ID: Q16539	
Isotype: IgG	Full Name: mitogen-activated protein kinase 14	
Immunogen Catalog Number: AG5115	Calculated MW: 360 aa, 41 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00437-1, p38 MAPK Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80821-5-PBS. Detection antibody: 80821-6-PBS. Standard: Ag5115. Range: 1.56-100 ng/mL.