

For Research Use Only

MCP-1/CCL2 Recombinant antibody, PBS Only (Detector)

Catalog Number:80847-6-PBS



Basic Information

Catalog Number: 80847-6-PBS	GenBank Accession Number: BC009716	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6347	CloneNo.: 241165F12
Source: Rabbit	UNIPROT ID: P13500	
Isotype: IgG	Full Name: chemokine (C-C motif) ligand 2	
	Calculated MW: 99 aa, 11 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00938-3, MCP-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80847-3-PBS. Detection antibody: 80847-6-PBS. Standard: Eg31375. Range: 0.313-20 ng/mL