

For Research Use Only

Ag11509 Recombinant antibody, PBS Only (Detector)

Catalog Number:80876-5-PBS



Basic Information

Catalog Number: 80876-5-PBS	GenBank Accession Number: BC050284	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 54915	CloneNo.: 241445G2
Source: Rabbit	UNIPROT ID: Q9BYJ9	
Isotype: IgG	Full Name: YTH domain family, member 1	
Immunogen Catalog Number: AG11509	Calculated MW: 559 aa, 61 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01146-2, YTHDF1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80876-4-PBS. Detection antibody: 80876-5-PBS. Standard: Ag11509. Range: 0.938-60 ng/mL.