

For Research Use Only

HCN2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 81025-4-PBS



Basic Information

Catalog Number:

81025-4-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001194

GeneID (NCBI):

610

UNIPROT ID:

Q9UL51

Full Name:

hyperpolarization activated cyclic
nucleotide-gated potassium channel
2

Calculated MW:

97 kDa

Purification Method:

Protein A purification

CloneNo.:

230398C8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00160-1, HCN2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81025-5-PBS. Detection antibody: 81025-4-PBS. Standard: Ag34038. Range: 0.313-40 ng/mL.