

For Research Use Only

TNFRSF14 Recombinant antibody, PBS Only (Detector)



Catalog Number: 81078-2-PBS

Basic Information

Catalog Number: 81078-2-PBS	GenBank Accession Number: BC002794	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 8764	CloneNo.: 240118B1
Source: Rabbit	UNIPROT ID: Q92956	
Isotype: IgG	Full Name: tumor necrosis factor receptor superfamily, member 14 (herpesvirus entry mediator)	
	Calculated MW: 30 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00364-1, TNFRSF 14 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81078-3-PBS. Detection antibody: 81078-2-PBS. Standard:null. Range: null.