

For Research Use Only

YAP1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 81087-1-PBS

Basic Information

Catalog Number: 81087-1-PBS	GenBank Accession Number: BC038235	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 10413	CloneNo.: 5H2
Source: Rabbit	UNIPROT ID: P46937	Affinity: $K_D=1.01 \times 10^{-10}$
Isotype: IgG	Full Name: Yes-associated protein 1, 65kDa	$K_{off}=4.24 \times 10^{-5}$
Immunogen Catalog Number: AG4510	Calculated MW: 504 aa, 54 kDa	$K_{on}=4.18 \times 10^5$

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00309-2, YAP1
Recombinant Matched Antibody Pair - PBS only.
Capture antibody: 81096-1-PBS. Detection
antibody: 81087-1-PBS. Standard: Ag4510. Range:
78.1-5000 pg/mL