

For Research Use Only

# YAP1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 81096-1-PBS

## Basic Information

<b>Catalog Number:</b> 81096-1-PBS	<b>GenBank Accession Number:</b> BC038235	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1mg/ml	<b>GeneID (NCBI):</b> 10413	<b>CloneNo.:</b> 1M10
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P46937	<b>Affinity:</b> $K_D=1.37 \times 10^{-10}$
<b>Isotype:</b> IgG	<b>Full Name:</b> Yes-associated protein 1, 65kDa	$K_{Off}=1.08 \times 10^{-4}$
<b>Immunogen Catalog Number:</b> AG4510	<b>Calculated MW:</b> 504 aa, 54 kDa	$K_{On}=7.87 \times 10^5$

## Applications

**Tested Applications:**  
Sandwich ELISA, Indirect ELISA

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP00309-1, YAP1 Recombinant Matched Antibody Pair - PBS only. Capture antibody: 81096-1-PBS. Detection antibody: 81099-1-PBS. Standard: Ag4510. Range: 78.1-5000 pg/mL



Sandwich ELISA standard curve of MP00309-2, YAP1 Recombinant Matched Antibody Pair - PBS only. Capture antibody: 81096-1-PBS. Detection antibody: 81087-1-PBS. Standard: Ag4510. Range: 78.1-5000 pg/mL