

For Research Use Only

PCNA Recombinant antibody, PBS Only (Capture)



Catalog Number: 81302-2-PBS

Basic Information

Catalog Number: 81302-2-PBS	GenBank Accession Number: BC000491	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5111	CloneNo.: 240306A6
Source: Rabbit	UNIPROT ID: P12004	
Isotype: IgG	Full Name: proliferating cell nuclear antigen	
Immunogen Catalog Number: AG0277	Calculated MW: 29 kDa/31 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00424-1, PCNA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81302-2-PBS. Detection antibody: 81302-4-PBS. Standard: Ag0277. Range: 0.78-100 ng/mL



Cytometric bead array standard curve of MP00424-3, PCNA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81302-2-PBS. Detection antibody: 81302-5-PBS. Standard: Ag0277. Range: 1.56-100 ng/mL