

For Research Use Only

P16-INK4A Recombinant antibody, PBS Only (Detector)

Catalog Number: 81373-3-PBS



Basic Information

Catalog Number: 81373-3-PBS	GenBank Accession Number: BC021998	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1029	CloneNo.: 241438D5
Source: Rabbit	UNIPROT ID: Q8N726	
Isotype: IgG	Full Name: cyclin-dependent kinase inhibitor 2A (melanoma, p16, inhibits CDK4)	
Immunogen Catalog Number: AG1328	Calculated MW: 16 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01135-2, P16-INK4A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81373-6-PBS. Detection antibody: 81373-3-PBS. Standard: Ag1328. Range: 0.781-100 ng/mL.